

## INSTRUMENT STERILIZER:-

### APPLICATIONS & KEY FEATURES

Brand	APS
Usage/Application	Hospital
Material	Stainless Steel
Material Grade	SS 304
Power Source	Electric
Voltage	110-220 v



### Product Description

**APS LAB INSTRUMENTS PVT LTD** are engaged in manufacturing and supplying the high quality assortment of **instrument sterilizer** under its brand name **APS** located in **delhi, india**.

instrument sterilizers electric (seamless), with hinged lid fitted with immersion rod and water outlet tap, seamless stainless steel construction. Fitted with automatic cutoff safety pin device, thermostat & tray lifting arrangement

**Note: Customized Solution Available On Request**