

LED OT LIGHT UNIQUE SERIES MODEL APS-OL-301

APS For fulfill our customer's requirements, we design a unique series with extra high luminance, brightness, high color purity and no radiation. This new model is fully designed for Large Modular Operation Theaters in single dome. It has all advance features of high quality surgical light.

Latest Modern technology extreme effective, excellent laminar flow, shadow less operation light with long life.

Features:

Light Color Pattern Adjustment:

The panel & remote is designed with graphical icons to allow user to quickly learn and easily operate:

- On/off finger touch switch
- Light intensity adjustment
- Light color pattern adjustment

Optional Camera & monitoring System

We can provide Camera (SD/HD) with Monitor also for modular operation theaters uses.

TECHNICAL DATA

| | |
|--------------------------|----------------|
| ➤ Luminance | 1,35,000Lux |
| ➤ Spot Diameter | 100 - 180 mm |
| ➤ Depth Illumination | 15 - 21 cm |
| ➤ Colour Temp. | 3200-5200 K |
| ➤ Colour Rendering Index | 91 Ra |
| ➤ Focusing | Manual |
| ➤ LED Average Life | >50000 Hrs. |
| ➤ Number of LED | 98 |
| ➤ Diameter of Light | 690 mm |
| ➤ Intensity Control | 10~100 Digital |
| ➤ Power Supply | 220V/50Hz AC |
| ➤ Power Consumption | 24 Watt |
| ➤ Special Function | Memory |



APS-OL-301

Technical Data

| | |
|------------------------|----------------|
| Intensity | 1,35,000 Lux |
| Spot Diameter | 100 - 180 mm |
| Depth Illumination | 15 - 21 cm |
| Colour Temp. | 4000-4850 K |
| Colour Rendering Index | 91 Ra |
| Focusing | Manual |
| LED Average Life | >50000 Hrs. |
| Number of LED | 98 |
| Diameter of Light | 690 mm |
| Intensity Control | 10~100 Digital |
| Power Supply | 220V/50Hz AC |
| Power Consumption | 55 Watt |
| Special Function | Memory |



CORPORATE OFFICE : C-9/36, BASEMENT, SECTOR-8, ROHINI, DELHI-110085

EMAIL ID : aps.inst2010@gmail.com, WEB: www.apsinstruments.com www.labequipments.net



AN ISO 9001 : 2008 & 14001 : 2004



MANUFACTURERS & EXPORTERS OF SCIENTIFIC INSTRUMENTS FOR LABS • HOSPITALS • SCHOOLS • COLLEGES & INDUSTRIES